INSTRUCTIONS FOR USE

ACTIVASORB® Depth Filter Sheets

Installation of ACTIVASORB® depth filter sheets

Before use, check the filter sheet conditions. Defective filter sheets should not be used. The sheets must be inserted with the unfiltered side (inlet/rough side) against the non-filtrating elements or frames and the filtrating side (outlet/fine side) against the filtrating elements of the filter.

NOTE

For the NOVOX® ST and OD plate and frame filter press, the filter sheets must be centered and cover all sides of the elements so that they protrude a few millimeters above the filter elements.

Wetting the filter sheets

- Gently press the plate and frame filter (approx. 50 % pressing power)
- Fill the filter with water (e. g. tap water or WFI) by slowly opening the inlet keeping the venting valves open
- Once purged, open the outlet valve to drain and close the venting valves, adjust the filter pressing power as required
- Apply back pressure of about 0.2 bar (2.9 psi) by closing the filter outlet slightly
- Rinse for 5 minutes to wet all filter sheets
- Stop the water supply and close the filter pack with full power

For manually operated, smaller filters, the wetting of the filter sheets can also be carried from the outside.

Rinsing the filter sheets

- Use 50 L/m² of a process-compatible solution (e. g. WFI) indicatively (might vary depending on ACTIVASORB® grade)
- The rinsing flow rate shall be approximately 350 L/m²×h

If rinsing is not desired or tolerated (insoluble products such as high percentage distillates or oils), the product to be filtered may be first recirculated to remove possible small fibers and/or diatomite particles.

Sterilization

- Slightly open (release) the plate and frame filter press (to ease the materials expansion due to the increasing temperature)
- Fill the system with water, slightly open all outlet valves
- Heat up the filter until the required sterilization temperature of 85 °C is reached at all outlet valves
- Maintain temperature for minimum 20 minutes
- Turn off the hot water supply
- Cool the filter down gradually to room temperature with cold sterile water while recirculating it into the system

To save energy, it is recommended in step k) and l) to heat up during recirculation.

NOTE

Sterilization is not fully completed until all parts to come into contact with the product have been sterilized. For alternative sterilisation methods, please contact your FILTRIX representative.

Filtration

- Slowly start filling the filter by opening the inlet valve and all venting valves (ensure complete venting), the outlet valve should be closed
- Once properly vented, open the outlet valve, close the venting valves and start the filtration (if required, the initial filtered volume may be separated from filtrate to avoid dilution with rinsing water)
- Adjust to the flow rate to 150 L/m²×h or accordingly to the process needs
- Ensure continuous and uniform flow rate
- Avoid any abrupt changes in filtration flow, pressure shocks, gas bubbles, as well as stop and go

NOTE

Observe the maximum recommended flow rate, operating pressure and temperatures. Further information can be found in the corresponding technical data sheet.

Regeneration

Once saturated and the filtration is no longer in accordance with the requirements, the depth filter sheets are to be considered exhausted and cannot be regenerated.

Replacement of filter sheets

The ACTIVASORB® depth filter sheets must be replaced as soon as the filtration-treatment requirements are no longer fulfilled. The filtered volumes vary depending on the products, process target, flow rate as well as ACTIVASORB® grade.

Note on disposal

The respective official regulations for disposal must be followed, depending on the filtered product. Uncontaminated sheets can be disposed of as non-hazardous waste. Please refer to your local environmental authority for more information.

Shelf life and storage

The filter sheets should be stored in their original packaging in a dry, odor-free and well-ventilated area. Temperature fluctuations can lead to condensation formation and should therefore be avoided. The filter sheets are intended for immediate use. Filter sheets for food & beverages should be used within 60 months from date of manufacture. Filter sheets for life science applications should be used within 36 months from date of manufacture. Further information can be found on the corresponding technical data sheet for filter media.

Technical Support

Contact your local FILTRIX representative or our experts at applications@filtrix.com. We will be happy to help you to find the best possible filtration solution for your product.

Remarks

All information is based on current knowledge and does not claim to be complete. No liabilities can be derived from this information. FILTRIX reserves the right to make changes in the course of technical improvement.

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