INSTRUCTIONS FOR USE

FIBRADISC® PLUS2 Depth Filter Modules

Installation of FIBRADISC® PLUS2 depth filter modules

Before use, check the filter modules conditions including gaskets. Defective filter modules should not be used. The modules must be inserted centrally according to their adapter and shall not get in contact with the housing at any point.

Rinsing the filter modules

- Use 50 L/m² of a process-compatible solution (e. g. WFI, ethanol)
- The rinsing flow rate shall be approximately 1.5 × L/m²×h to the given filtration flow for the corresponding product to be filtered
- Rinse in the filtration direction (forward flow)

If the rinsing is not desired or tolerated (insoluble products such as high percentage distillates or oils), the product to be filtered may be first recirculated to remove possible small fibers and/or diatomite particles.

Sterilization

- Fill the system with water, slightly open all outlet valves and increase the temperature gradually
- Heat up the filter until the required sterilization temperature of 85 °C is reached at all outlet valves
- Maintain temperature for minimum 20 minutes
- Turn off the hot water supply
- Cool down gradually the filter to room temperature with cold sterile water while recirculating it into the system

To save energy, it is recommended in d) and f) to heat up during circulation.

NOTE

Sterilization is not fully completed until all parts that come into contact with the product have been sterilized. For alternative sterilisation methods, please contact your FILTROX representative.

Filtration

- Slowly start filling the filter by opening the inlet valve and all venting valves (ensure complete venting), the outlet valve should be closed
- In case of beer or carbonated products, keep enough back pressure regulating the venting valve to avoiding foam formation during filling
- Once properly vented, open the outlet valve, close the venting valves and start the filtration (if required, the initial filtered volume may be separated from the filtrate to avoid dilution with rinsing water)
- Adjust the flow rate to the installed filter module grade, process requirements and filtration targets for the product to be filtered
- Ensure a continuous and uniform flow rate
- Avoid any abrupt changes in the filtration flow, pressure shocks, gas bubbles as well as stop and go

NOTE

Observe the maximum recommended flow rate, operating pressure and temperatures. Further information can be found in the corresponding technical data sheet.

Back flushing

As soon as differential pressure exceeds the starting differential pressure approx. 0.3 – 0.5 bar, a back flush can be started to help extend the module life and yield.

- Fill the filter with cold water or push out the filtered product with water
- Rinse forward for 10 – 15 min. with cold water at 1.5 × the given nominal flow rate (if differential pressure allows it)
- Rinse back flush for 5 – 10 min. at 0.4 bar with warm water at 40 °C (max. 0.5 bar / 7 psi at 40 °C)
- Rinse again forward for 5 min. with cold water and check the differential pressure

After a successful back flush, the differential pressure should be lower than mentioned above. If so, proceed with step p) and then restart the filtration. If the pressure has not decreased, continue backflushing.

Regeneration

- Fill the filter with cold water
- Rinse for 10 – 15 min. with cold water at 1.5 × the given nominal flow rate (if differential pressure allows it)
- Rinse for 10 – 15 min. with warm water (max. 45 °C) at 1.5 × the given nominal flow rate (if differential pressure allows it)
- If a sterilization is required, follow the sterilization process in e)
- Cool the filter down gradually to room temperature with cold sterile water while recirculating it into the system

NOTE

Do not exceed the max. permissible differential pressure anytime.

Filter modules wet storage (in the housing)

After regeneration, <1 week in a 2% citric acid solution. >1 week in a 0.2% SO₂ or peracetic acid solution approved for food and beverages applications. Replace the solution every week. Before use, the modules must be rinsed to ensure chemical elimination and sterilized.

Replacement of filter modules

The depth filter modules must be replaced as soon as the maximum permissible differential pressure has been reached. Please refer to the specified pressure limits in the corresponding data sheet.

Note on disposal

The respective official regulations for disposal must be followed depending on the filtered product. Uncontaminated modules can be disposed of as non-hazardous waste. Please refer to your local environmental authority for more information.

INSTRUCTIONS FOR USE

FIBRADISC® PLUS2 Depth Filter Modules

Shelf life and storage

The filter modules should be stored in their original packaging in a dry, odor-free and well-ventilated area. Temperature fluctuations can lead to condensation formation and should therefore be avoided. The filter modules are intended for immediate use. Filter modules for food & beverages should be used within 60 months from date of manufacture. Filter modules for life science applications should be used within 36 months from date of manufacture. Further information is given in the corresponding technical data sheet for filter media.

Technical Support

Contact your local FILTRIX representative or our experts at applications@filtrix.com. We will be happy to help you to find the best possible filtration solution for your product.

Remarks

All information is based on current knowledge and does not claim to be complete. No liabilities can be derived from this information. FILTRIX reserves the right to make changes in the course of technical improvement.

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